Kuglen Iris Hook K-shaped tip

Ref	Tip	
554000	Angled	
554001	Straight	



Jaffe-Knolle Iris Hook Iris Hook

Ref	Tip
554002	Angled
554003	Straight



Fukasaku Small Pupil Snapper Hook With smoothly polished notch and sharp inner edge

Ref Tip 554004 Angled



Lens Manipulator Fenzl hook

Ref	Tip
554005	Straight
554006	Angled



Lens Manipulator Y nucleus rotator angled


Ref	Tip	
554007	Vertical	
554087	Horizontal	



Lens Manipulator V-shaped tip

Ref	Tip
554008	Angled
554009	Straight





Lens Manipulator

1.8 mm diameter - Angled

Ref	Overall length
554010	120 mm





Lens Manipulator

1.1 m diameter - Curved

Ref	Overall length
554011	120 mm





Nucleus Manipulator Drysdale tip

Ref	Overall length
554012	120 mm





Manipulator Clayman Spatula, Flat

Ref	Overall length
654119	120 mm





Manipulator Sinskey, tip 0,2mm

Ref	Altitude
654120	0.7 mm
654121	1.1 mm





Lieberman Microfinger

Ref	Overall length
654122	115 mm
100	
,	



Manipulator Cylinder

Ref	Overall length
654112	120 mm



Manipulator govan

Ref	Overall length
654113	125 mm



Nucleus Manipulator Mushroom manipulator for iris retraction

Ref	Overall length
554013	120 mm



Nucleus Manipulator Connor wand

Ref	Overall length	
554014	120 mm	





Nucleus Manipulator Connor wand

Ref	Overall length
554015	120 mm
CONTRACTOR OF THE PARTY OF THE	



Nucleus Manipulator Connor wand

Ref	Overall length	
554016	120 mm	
		-



Iris Hook H-shaped tip

Ref	Tip
554017	Straight
554018	Angled



Kuglen Iris Hook "S" shaped 0.5 mm pointed hook

Ref	Overall length
554066	120 mm



Kuglen Iris Hook Curved 0.5 mm hook

Ref	Overall length
554067	120 mm



Lund Iris Hook 14 mm tip, 45° Angled

Ref	Overall length
654088	120 mm



Manipulator Lester Angled


Ref	Overall length	
654097	120 mm	



Manipulator Sinskey 0,2mm hooked angle

Ref	Overall length
654100	120 mm





Phaco Chopper Angled shaft, 1.8 mm tip

Ref	Overall length
554062	Left
554063	Right



Phaco Chopper Chopper and manipulator in one

Ref	Hand	
554064	120 mm	



Phaco Chopper Special inner blade is offset 45°

Ref	Overall length
554065	120 mm




Ref	Overall length
654114	115 mm

Koch Stop & Chop Manipulator

Ref	Overall length
554085	115 mm

Repositor Semi-sharp tip 0.7 mm



Repositor Semi-sharp tip 0.5 mm

Ref	Overall length
554028	120 mm
_//	



Ref	Overall length
554071	120 mm



ICL Manipulator

554073	Angled
554074	Curved



Repositor Doubled-ended iris spatula




Lebuisson Manipulator

Ref	Overall length
654107	120 mm

Lasik Flap lifter/sinskey

Ref	Overall length
654124	115 mm





LASIK Flap Separator/Elevator

Ref	Overall length
554503	130 mm



Lasik Flap Lifter/Elevator Binder

Ref	Overall length
554504	115 mm





Lenticule Separator/Elevator

Ref	Overall length
654108	130 mm
P	





Double-ended Manipulator Drysdale/H-shaped tip

Ref	Overall length
654085	130 mm



Double-ended Manipulator Drysdale/Y pusher

Ref	Overall length
654086	125 mm



Double-ended Manipulator Cylinder/Y pusher

Ref	Overall length
654087	125 mm



Double-ended Manipulator Drysdale/ Fukasaku small pupill snapper

Ref	Overall length
654089	125 mm

MAA.07T ⊕ 654000 C €

Double-ended Manipulator Agarwal Phaco Chopper/Fukasaku small pupil snapper

Ref	Overall length
654090	125 mm



Double-ended Manipulator Sinskey/H-shaped tip

Ref	Overall length
654091	125 mm



Double-ended Manipulator Drysdale/ Koch Chopper

Ref	Overall length
654096	125 mm



Double-ended Manipulator Fukasaku Small Pupil Snapper/cylinder

Ref	Overall length
654104	125 mm



Double-ended Manipulator Koch Stop & Chop/Phaco Chopper

Ref	Overall length
654106	125 mm



Double-ended Manipulator Sinskey/Lester Hook

Ref	Overall length
654109	130 mm

Double-ended Manipulator Lund Iris Hook/Phaco Chopper

Ref	Overall length
654110	130 mm



Double-ended Manipulator Mushroom/Phaco Chopper

Ref	Overall length
654111	130 mm



Double-ended Manipulator Sweep/Y-horizontal

Ref	Overall length
654118	130 mm



Double ended reversed and sinskey

Ref	Overall length
654123	125 mm



Double ended drysdale and chopper

Ref	Overall length
554501	125 mm



Fukasaku & Sinskey

Ref	Overall length
554502	125mm



Buratto Flap Protector Ideal for retracting the flap

Ref	Overall length
554041	120 mm



Spatula Corneal epithelium spatula

Ref	Overall length
554082	120 mm



Epi-peeler

For manipulating the epithelium in LASIK procedures

Ref	Overall length
554083	120 mm



Hoe Type Epi-peeler Square shaped

Ref	Overall length
554084	120 mm





Helveston Retractors

Ref	Tip length
554042	9 mm
554043	7 mm
554044	11 mm



Kansas Nucleus Vectis Double open ends

Ref	Overall length
554051	130 mm



Kansas Nucleus Vectis Double solid ends

Ref	Overall length
554052	130 mm



Kansas Nucleus Bisector For manual small incision surgery

Ref	Overall length
554053	120 mm



Lens Loop

Ref	Loop
554034	5x7 mm
554035	3x8 mm

