

MARKERS & GAUGES

Markers - Calipers - Gauges

Lasek Alcohol Well

Ref	Diameter
554036	7.0 mm
554037	8.0 mm
554038	9.0 mm
654039	9.5 mm



Optical Zone - Target Centre 4.0 mm diameter

Ref	Overall length
558000	111 mm



CK Lasik Marker - Target Centre 10 mm diameter

Ref	Overall length
558001	119 mm



Mendez Degree Gauge - Calibrated every 10 degrees from 0° to 180° External diameter 19.5 mm Internal diameter 12.5 mm

Ref	Overall length
558002	112 mm



Lasik Corneal Marker Two blunt fins and separated at 45° angle

Ref	Overall length
558004	119 mm



MARKERS & GAUGES

Markers - Calipers - Gauges

Tsuneoka micro incision gauges

Set of 8 stainless steel blades

Ranges from 1.8 mm to 2.5 mm in 1 mm increments

Ref	Overall length
558007	50 mm



Castroviejo Caliper

Scale reads on one side

Measure from 0 mm to 20 mm in 0.5 mm increments

Ref	Overall length
558009	88 mm




Braunstein Fixed Caliper

Stainless steel blades

Size from 1.8 mm to 2.5 mm in 1 mm increments

Ref	Overall length
558008	70 mm



Scleral Marker - Double ended

Aphakic 3.5/4.0 mm

Ref	Overall length
558003	115 mm



Braunstein Fixed Caliper - Fine tips mark sclerotomy sites for vitreous surgery

Aphakic 3.5/4.0 mm

Ref	Overall length
558005	68 mm



MARKERS & GAUGES

Markers - Calipers - Gauges

Braunstein Fixed Caliper - Fine tips mark sclerotomy sites for vitreous surgery Aphakic 3.5/4.0 mm

Ref	Overall length
558006	69 mm



Braunstein Fixed Caliper Straight *

Ref	Overall length
658011	75 mm

* only available in IVT packs



Fixation Marker 3.5 *

Ref	Overall length
658013	70 mm

* only available in IVT packs



IVT Marker Small *

Ref	Overall length
658012	110 mm

* only available in IVT packs



Plastic Braunstein Fixed Caliper *

Ref	Overall length
558010	80 mm

* only available in IVT packs



MARKERS & GAUGES

Markers - Calipers - Gauges

Braunstein Fixed Caliper *

Fine tips mark sclerotomy sites for vitreous surgery

Ref	Overall length
658005	90 mm

* only available in IVT packs

