

MICRO INSTRUMENTS

Micro Instruments

Grasping Forceps smooth, 45°

Ref	Edge
556035	without
556037	with



Micro Forceps Straight forceps, with 1x2 teeth

Ref	Gauge
556031	20 g
556032	23 g
656032	25 g



IOL Cutter 19G

Ref	Overall length
653057	114 mm



Scissors, Curved shaft sharp tips

Ref	Gauge
556017	23 g
656017	25 g



Scissors, Curved long shaft sharp tips

Ref	Gauge
556018	23 g
656018	25 g



MICRO INSTRUMENTS

Micro Instruments

Trabeculectomy Punch

Ref	Tip length
554020	0.75 mm
554021	1.0 mm
554022	1.5 mm
554023	2.0 mm

