

SA60-2200-10 Fig. 1

SA60-2200-20 Fig. 2

SA60-2200-30 Fig. 3

WILDER (NETTLESHIP)

4 mm Ø

Tränenweg-Dilatator

Lacrimal dilator

Dilatateur pour voies lacrimales

Dilatador lacrimal



SA60-2201-10 Fig. 1

SA60-2201-20 Fig. 2

SA60-2201-30 Fig. 3

WILDER


3 mm Ø

Tränenweg-Dilatator

Lacrimal dilator

Dilatateur pour voies lacrimales

Dilatador lacrimal



SA60-2202-00

CASTROVIEJO

Dilatator, doppelendig

Lacrimal dilator, double ended

Dilatateur pour voies lacrimales

double

Dilatador lacrimal doble



SA60-2203-00

MULDOON

Dilatator

Lacrimal dilator

Dilatateur pour voies lacrimales

Dilatador lacrimal



SA60-2206-00

RUEDEMANN

Dilatator für Säuglinge

Lacrimal dilator for infants

Dilatateur pour nourrissons

Dilatador lacrimal para lactantes



SA60-2207-00

MODELL GIESSEN

Tränenwegsonde, abgewinkelt

Lacrimal probe, angled

Sonde a voies lacrymales, coude

Estilete lacrimal, angulada



GEFÄSSKLEMMEN

SERREFINES

SERRE-FINES

SERRA-FINAS

SA60-2260-35 JOHN HOPKINS

3,8 cm

Gefäßklemme, gerade
Serrefine, straight Serre-
fine, droit
Serra-fina, recta



SA60-2261-35 JOHN HOPKINS

3,8 cm

Gefäßklemme, gebogen
Serrefine, curved
Serre-fine, courbe
Serra-fina, curva



KLEMMEN

CLAMPS


PINCES-CLAMPS

PINZAS-CLAMPS

SA60-2265-95 BACKHAUS

9,5 cm

Tuchklemme
Towel clamp
Pince a etoffe
Pinza de campo



SA60-2267-50 JONES

Tuchklemme
Towel clamp, cross action
Pince a compreses
Pinza para compresas



SA60-2272-10
HARTMANN-MOSQUITO


10 cm

Arterienklemme, gerade
Mosquito forceps, straight
Pince hemostatique, droite
Pinza hemostatica, recta

SA60-2273-10
HARTMANN-MOSQUITO

10 cm

Arterienklemme, gebogen
Mosquito forceps, curved
Pince hemostatique, courbe
Pinza hemostatica, curva



SA60-2274-12
HALSTED-MOSQUITO


12,5 cm

Arterienklemme, gerade
Mosquito forceps, straight
Pince hemostatique, droite
Pinza hemostatica, recta

SA60-2275-12
HALSTED-MOSQUITO

12,5 cm

Arterienklemme, gebogen
Mosquito forceps, curved
Pince hemostatique, courbe
Pinza hemostatica, curva





SA60-2285-15

GASS

1,5 mm

Sklerotomie Stanze
Sclerotomy punch
Emporte-pièce
Punzón





SA60-2286-10

KELLY

1,0 mm

Descement Membran Stanze
Descement's membrane punch
Emporte-pièce
Punzón




SA60-2288-10 1,0 mm

SA60-2288-15 1,5 mm

HOLTH

Kapselstanze
Corneoscleral punch
Emporte-pièce capsulaire
Punzón capsular





SA60-2290-10

LUNTZ-DODICK

1,0 mm



Trabekulektomie Stanze
Trabeculectomy punch
Punch trabeculectomie
Punzón para trabeculectomía




STANZEN
PUNCHES

EMPORTE-PIÈCES
PUNZONES

SA60-2293-30 3,0 mm
SA60-2294-40 4,0 mm
KERRISON
Knochenstanze
Rongeur
Emporte-pièce
Pinza-gubia



2SA60-296-20 2,0 mm
SA60-2296-30 3,0 mm
CITELLI
Knochenstanze
Rongeur
Emporte-pièce
Pinza-gubia



SA60-2298-20
BEYER
2,0 mm
Knochenstanze
Rongeur
Emporte-pièce
Pinza-gubia

