


Normale Größe – Regular size – Taille normal – Tamaño normal

SA60-4040-41

VANNAS

8,5 cm

gerade, 7 mm Blattlänge
straight, 7 mm blades
droit, lames 7 mm
recto, hojas 7 mm de largo




SA60-4041-41

VANNAS

8,5 cm

gebogen, Blattlänge
7 mm curved, 7 mm blades
courbe, lames 7 mm
curvo, hojas 7 mm de largo




SA60-4042-41

VANNAS

8,5 cm

leicht aufgewinkelt
slightly angled on flat
légèrement angulaire
ligeramente angulado



SA60-4043-41

VANNAS

8,5 cm

stark aufgewinkelt
strong angled on flat
très angulaire
muy angulado




SA60-4044-41

VANNAS

8,0 cm

kniegebogen
angled to side
coudé
angulado




Extra fein – Extremely delicate – Extra fine – Delicadísima

SA60-4046-41
VANNAS

8,5 cm


gerade, 5 mm Blattlänge
straight, 5 mm blades
droit, lames 5 mm
recto, hojas 5 mm de largo



SA60-4047-41
VANNAS

8,5 cm


gebogen, Blattlänge 5 mm
curved, 5 mm blades
courbe, lames 5 mm
curvo, hojas 5 mm de largo



SA60-4048-41
VANNAS

8 cm


kniegebogen
angled to side
coudé
angulado



SA60-4049-41
VANNAS

8,5 cm

leicht aufgewinkelt
slightly angled on flat
légèrement angulaire
ligeramente angulado



VANNAS-ARTIGE SCHEREN

VANNAS TYPE SCISSORS

CISEAUX VANNAS

TIJERAS TIPO VANNAS

SA60-4050-41

MODELL TUEBINGEN

9 cm

gerade, Blattlänge 8 mm
straight, 8 mm blades
droit, lames 8 mm
recto, hojas 8 mm de largo



SA60-4051-41

MODELL TUEBINGEN

9 cm

gebogen, Blattlänge 8 mm
curved, 8 mm blades
courbe, lames 8 mm
curvo, hojas 8 mm de largo



SA60-4053-41

MODELL TUEBINGEN

8,5 cm

stark gebogen, Blattlänge 9 mm
strong curved, 9 mm blades
très courbe, lames 9 mm
muy curvo, hojas 9 mm de largo




SA60-4052-41

GILL

8,5 cm

aufgewinkelt
angled on flat
angulaire
angulado



SA60-6281-41

GILLS-WELSH

8,5 cm

aufgewinkelt
angled on flat
angulaire
angulado




SA60-4060-61

CASTROVIEJO-VANNAS

12 cm

gerade, Blattlänge 6 mm
straight, 6 mm blades
droit, lames 6 mm
recto, hojas 6 mm de largo




SA60-4061-61

CASTROVIEJO-VANNAS

12 cm


gebogen, Blattlänge 6 mm
curved, 6 mm blades
courbe, lames 6 mm
curvo, hojas 6 mm de largo



SA60-9835-61


VANNAS

12,5 cm
gebogen
curved
courbe
curvo




SA60-4065-43

8,5 cm
gebogen, leicht stumpf, Blattlänge 6 mm
curved, slightly blunt, 6 mm blades
courbe, mousse, lames 6 mm
curvo, romo, hojas 6 mm de largo




SA60-4067-43

9 cm
gebogen, leicht stumpf, Blattlänge 6,5 mm
curved, slightly blunt, 6,5 mm blades
courbe, mousse, lames 6,5 mm
curvo, romo, hojas 6,5 mm de largo



SA60-4070-41

8,5 cm
Trabekular-Schere, gerade,
Blattlänge 6,5 mm
Trabecular scissors, straight,
6,5 mm blades
Ciseaux à trabéculaires, droit,
lames 6,5 mm
Tijera trabeculares, recta,
hojas 6,5 mm de largo



SA60-4071-41

Trabekular-Schere, gebogen,
Blattlänge 6,5 mm
Trabecular scissors, curved,
6,5 mm blades
Ciseaux à trabéculaires, courbe,
lames 6,5 mm
Tijera trabeculares, curva,
hojas 6,5 mm de largo

