




Optionaler Thermolock 410.47.01.0

Deckel M

EasyStop 431.03.01.X	Permanentfilter 433.03.01.X	Papierfilter 432.03.01.X	Entsorgungsdeckel 437.04.01.X
-------------------------	--------------------------------	-----------------------------	----------------------------------

rot 1 silber 2 blau 3 grün 4 gold 5 schwarz 6 grau 7

Wanne

Größe	Außenabmessung	Innenabmessung [max. Beladehöhe]	Artikelnummer
M 100	290 x 290 x 100	253 x 269 x [73]	430.05.10.2.XX.00
M 135	290 x 290 x 135	253 x 269 x [108]	430.05.13.2.XX.00
M 150	290 x 290 x 150	253 x 269 x [123]	430.05.15.2.XX.00
M 170	290 x 290 x 170	253 x 269 x [143]	430.05.17.2.XX.00
M 220	290 x 290 x 220	253 x 269 x [193]	430.05.22.2.XX.00
M 270	290 x 290 x 270	253 x 269 x [243]	430.05.27.2.XX.00



ohne Thermolock 06



ohne Thermolock 10