


Just One Use // The acronym « *Just One Use* » affirms the *single use characteristic of our instruments*. **Safety** : reduce the risk of cross-contamination. **Economy** : eliminate the costs of repair and sterilization. **Performance** : manufactured from medical grade stainless steel, our single use instruments have the same grip, weight and balance of reusable instruments. **Convenience** : reduce theatre preparation time, presentation in sets or individually.

CAPSULORHEXIS / capsulorhexis / kapsulorhexis




SC33M

MICS CAPSULORHEXIS FORCEPS CROSS ACTION

Pince Capsulorhexis Action Croisée, Courbe
MICS Kapsulorhexispinzette Cross Action (x-Action)


Length : 115 mm
Handle length : 84 mm
50° angled
1.8 mm incision



CAPSULORHEXIS ROUND BODIED CURVED / Pince Capsulorhexis, Manche Rond, Courbe / Kapsulorhexispinzette, gebogen, runder Griff

Length : 105 mm
Handle length : 77 mm
45° angled platform with notched tips 0.3 mm
1.8 mm incision

SC34M

CAPSULORHEXIS

MICS CAPSULORHEXIS FORCEPS CURVED / Pince Capsulorhexis MICS, Courbe / MICS Kapsulorhexispinzette, gebogen



Length : 105 mm
Handle length : 77 mm
45° angle
12 mm platform
0.3 mm tips
1.8 mm incision

SC29M

CAPSULORHEXIS

MICS CAPSULORHEXIS FORCEPS ROUND BODIED CURVED / Pince Capsulorhexis, Manche rond / Kapsulorhexispinzette, runder Griff

Length : 105 mm
Handle length : 77 mm
50° angle 0.2 mm tip
12 mm platform
1.8 mm incision

SC134M

CAPSULORHEXIS

MICS CAPSULORHEXIS FORCEPS STRAIGHT / Pince Capsulorhexis MICS Droite, Kystitome / MICS Kapsulorhexispinzette, gerade

Length : 105 mm
Handle length : 76 mm
45° angled 0.3 mm notched tips, with 11 mm platform
1.8 mm incision

SC133M

CAPSULORHEXIS
