

Just One Use // Das Akronym « Just One Use » bestätigt die einzigartige Verwendung unserer Instrumente. Sicherheit :


Reduziert das Risiko einer Kreuzkontamination. **Wirtschaftlichkeit :** Kosten für Reparatur und Sterilisation entfallen.

Leistung : Verwendung von medizinischem Stahl, Instrumente mit einzigartiger Ergonomie, Gewicht und Beschaffenheit, die äquivalent zu sterilisierbaren Instrumenten ist. **Komfort :** Verkürzung der Vorbereitungszeit, erhältlich als Sets oder Einzelinstrumente.

MANIPULATORS & CHOPPERS

Manipulateurs & Choppers


Manipulatoren & Chopper



CHOPPER AND DIALLER / Double Manipulateur Chopper/Y / Chopper und Spalter

Length : 120 mm
Handle length : 76 mm
45° angled ends
One rounded tip and one 1 mm Y shaped end
Round bodied
Double ended


SC39



DENIS CLEAVER / Cliver de Denis / Denis Spalter

Handle length : 95 mm
0,3mm Curved Cleaver
Tip angled left 30° from 11mm 45° angled shaft
Round bodied


SC97



GREENS CHOPPER / Chopper de Greens / Greens Chopper

Length : 115 mm
Handle length : 95 mm
30° angled top curved rounded hook
Round bodied

SC45



GOVAN MANIPULATOR / Manipulateur Govan / Govan Manipulator

Length : 100 mm
Handle length : 75 mm
0.4 mm x 0.2 mm x 10.4 mm flattened tip
60° angle pointed tip
Round bodied


SC100



GOVAN MANIPULATOR AND PHACO CHOPPER / Chopper Phaco & Manipulateur Govan / Govan Manipulator & Phako Chopper

Length : 120 mm
Handle length : 78 mm
60° pointed end and 90° blunt end
Round bodied
Double ended

SC102



IRIS REPOSITOR / Manipulateur à iris / Iris Repositor

Length : 130 mm
Handle length : 98 mm
Round bodied


SC63



JONES CHOPPER / Chopper Jones / Jones Chopper

Length : 120 mm
Handle length : 95 mm
45° angled pointed tip
Round bodied

SC90



KOCH CHOPPER / Chopper Koch / Koch Chopper

Length : 120 mm
Handle length : 95 mm
45° angled 1 mm rounded tip
Round bodied


SC44



KOCH PHACO SPATULA / Spatule de Koch / Phako Spatel nach Koch

Length : 120 mm
Handle length : 96 mm
Slight curved spatula
Round bodied

SC143



KUGLEN HOOK / Crochet de Kuglen / Häkchen nach Kuglen

Length : 115 mm
Handle length : 95 mm
45° angled blunt, flat hooked tip
Round bodied

SG15



LEBUISSON MANIPULATOR / Manipulateur Lebuisson / Lebuisson Manipulator

0.40 mm x 0.50 mm x 1.00 mm spatula tip
Ø0.3 mm curved 12.6 mm
90° angled shaft
Length : 115 mm
Handle length : 90 mm
Round bodied

SC106



LENS DIALER (SINSKEY HOOK) / Crochet de Sinskey / Linsenhäkchen nach Sinskey

Length : 120 mm
Handle length : 95 mm
35° angled
Round bodied


SC46



LESTER HOOK / Crochet de Lester / Lester Haken

Handle length : 95 mm
Ø0.25 mm ball on 0.4 mm finger
10mm 45° angled shaft
Length : 115 mm
Handle length : 100 mm
Round bodied


SC146



PHACO CHOPPER AND MUSHROOM MANIPULATOR DOUBLE ENDED / Manipulateur Double Chopper-Champignon / Phaco Chopper und Mushroom Manipulator (Kombi-Instrument)

Length : 120 mm
Handle length : 70 mm
45° chopper
Round bodied
Double ended


SC147



PHACO CHOPPER / Chopper Phaco / Phako Chopper

Length : 120 mm
Handle length : 95 mm
40° angled round end
Round bodied

SC43



PHACO CHOPPER CLAW / Chopper Claw, fin / Phako Claw Chopper (krallen förmig, fein)

Length : 120 mm
Handle length : 95 mm
40° angled tip, hooked tip
Round bodied

SC42




PISACANO ROTATOR AND NEUHANN MANIPULATOR DOUBLE END / Manipulateur double Pisacano-Neuhann / Pisacano-Neuhann Manipulator (Kombi-instrument)

Length : 100 mm
Handle length : 75 mm
60° spatula end and 30° round end
Round bodied
Double ended


SC101

SINSKEY HOOK AND GOVAN MANIPULATOR DOUBLE ENDED / Manipulateur double Sinskey-Govan / Sinskey Håkchen und Govan Manipulator (Kombi-Instrument)

Length : 120 mm
Handle length : 76 mm
60° manipulator and 45° hooked end
Round bodied
Double ended

SC144

MUSHROOM NUCLEUS MANIPULATOR / Manipulateur champignon / Mushroom Nukleus Manipulator

Length : 120 mm
Handle length : 95 mm
45° angled mushroom shaped tip
Round bodied


SC48

Y SHAPED NUCLEUS ROTATOR / Manipulateur de Noyau Y, Vertical / Y-Nukleus Rotator

Length : 120 mm
Handle length : 95 mm
40° angled tip Y shaped end
Round bodied

SC41



CHOPPER AND KUGLEN HORIZONTAL DOUBLE ENDED / Manipulateur double Chopper/Kuglen / Kombi Manipulator Chopper/Kuglen

Length : 127 mm
Handle length : 70 mm
60° spatula shaped chopper
40° horizontal kuglen hook



NAGAHARA CLEAVER / Chopper Nagahara / Nagahara Chopper

Length : 117 mm
Handle length : 91 mm



DRYSDALE NUCLEUS SPATULA / Manipulateur Drysdale / Drysdale Nukleus Manipulator

Length : 110 mm
Handle length : 95 mm
45° angled paddle shaped tip
Round bodied

