

Just Quality // Quality of manufacture and irreproachable ethics.

Our instruments are CE marked and manufactured in compliance with current international standards. All our products are subject to rigorous quality testing throughout the design and manufacturing. We inspect each surgical instruments that we manufacture, and every instrument is hand finished. We strive to operate and behave in an ethical and responsible way.

SPECULUMS / Blépharostats / Lidsperrer



 **WIRE | MECHANICAL | PAEDIATRIC**

KRATZ CLOSED LOOP SPECULUM

Blépharostat temporal de Barraquer

Lidsperrer nach Barraquer, feiner gebogener Draht temporal offen

Length : 32 mm
50° angle at base



BROAD GAUGE WIRE SPECULUM / Blépharostat fils large / Breiter Gauge Draht Lidsperrer


Length : 50 mm

Width : 40 mm

SC01



WIRE



CLOSED LOOP WIDE SPECULUM / Blépharostat anses fils large / Spekulum/Lidsperrer mit geschlossenen Valven


Length : 44 mm

Width : 42 mm

SC07



WIRE




HEAVY GAUGE WIRE SPECULUM OPEN LOOP / Blépharostat fil, grand diamètre / Drahtlidsperrer groß (Heavy Gauge), offen

Length : 72 mm

Overall width : 58 mm

SC124



WIRE




KRATZ BARRAQUER CLOSED LOOP / Blépharostat de Barraquer, fenêtré / Lidsperrer nach Kratz , gefenstert

Length : 40 mm

Width : 42 mm

SC09




WIRE




KRATZ BARRAQUER CLOSED LOOP SPECULUM / Blépharostat Kratz-Barraquer, fenêtré / Lidsperrer Kratz-Barraquer, gefenstert

Length : 43 mm
50° angle at base

SC132




WIRE


KRATZ BARRAQUER WIRE SPECULUM OPEN LOOP / Blépharostat Kratz-Barraquer, anses fil / Kratz-Barraquer Lidsperrer, feiner gebogener Draht temporal offen

Length : 47 mm
Width : 31 mm
50° angle at base

SC11



WIRE



KRATZ OPEN LOOP WIRE SPECULUM (BLUE) / Blepharostat de Kratz (bleu) / Kratz Lidsperrer mit offenen Valven (Blau)

Length : 44 mm
Overall width : 30 mm
50° angle at base

SC107



WIRE



WIRE SPECULUM OPEN LOOP / Blepharostat Kratz-Barraquer fil, ouvert / Kratz-Barraquer Spekulum mit offenen Valven

Length : 38 mm
Overall width : 39 mm
Wire width : 0.9 mm

SC150



WIRE




KRATZ SOLID BLADED WIRE SPECULUM / Blépharostat Kratz-Barraquer, anses pleines (Chrome) / Kratz-Barraquer Lidsperrer, Valven aus Chrom (solid)

Length : 42 mm
Width : 30 mm
65° angle at base

SC05




WIRE




TEMPORAL SPECULUM / Blépharostat temporal de Kratz / Lidsperrer nach Kratz, feiner gebogener Draht temporal offen

Length : 50 mm
Width : 42 mm

SC12



WIRE




BANGERTER ADULT SPECULUM / Blepharostat de Bangerter / Lidsperrer für Erwachsene nach Bangerter

Length : 50 mm
sliding screw adjustable
Blade length : 15 mm

SC03



MECHANICAL




LIEBERMAN ADJUSTABLE SPECULUM, STAINLESS STEEL / Blépharostat Lieberman mécanique, anses fil / Mechanisches verstellbarer Lidsperrer nach Lieberman

Overall Length : 81 mm
Blade Length : (R) 14.5 mm / (L) 14.2 mm
Radius of blades : 2.5 mm
Thickness of wire : 0.9 mm

SC121



MECHANICAL




BLADED ADJUSTABLE SPECULUM / Blépharostat mécanique, anses pleines / Lidsperrer nach Barraquer, Valven aus Vollmaterial

Length : 75 mm
Blade length : 15 mm

SR10



MECHANICAL



CLARKE SPECULUM / Blépharostat de Clarke / Clarke Lidsperrer

Length : 80 mm
Open loop : 12 mm

SC126



MECHANICAL



KRATZ SCREW ADJUSTABLE WIRE SPECULUM OPEN LOOP / Blépharostat mécanique Kratz, anses fil / mechanisches verstellbarer Lidsperrer nach Kratz (feiner gebogener Draht) temporal offen

Length : 78 mm
Blade length : 15 mm

SR11



MECHANICAL



SOLID BLADED SCREW ADJUSTABLE SPECULUM / Blépharostat mécanique / Mechanisch feststellbarer Lidsperrer mit Valven aus Vollmaterial

Length : 85 mm
Blade length : 14 mm

SC08



MECHANICAL




KRATZ PAEDIATRIC CLOSED LOOP 10 MM SPECULUM / Blépharostat pédiatrique de Kratz / Paediatrischer Lidsperrerr nach Kratz 10mm

Length : 29 mm
Width : 30 mm
10 mm closed loop

SC10



PAEDIATRIC



NEONATE PAEDIATRIC SPECULUM / Blépharostat néonatal Necker / Neonatal paediatrischer Lidsperrerr

Length : 25 mm
Width : 20 mm
4 mm closed loop
80° angle at base

SC103



PAEDIATRIC




NEONATE 4 MM SPECULUM / Blépharostat néonatal 4 mm / Neonatal Lidsperrerr 4 mm

Length : 25 mm
Width : 31 mm
20-30 mm adjustable
3 mm blade

SC15



PAEDIATRIC



PAEDIATRIC 6 MM SPECULUM / Blépharostat pédiatrique 6mm / Paediatrischer Lidsperrerr 6mm

Length : 25 mm
Width : 35 mm
6 mm blade

SC14



PAEDIATRIC

