

Just Innovation // Continuous improvement and innovation are our driving force.

Our research and development team is constantly monitoring the market in order to offer surgical instruments in line with expressed needs. Our know-how is at the disposal of the users for any request for development of new or custom made instruments.

CURETTES & CLAMPS

Curettes & Pincers


Küretten & Klemmen



CHALAZION CLAMP MEDIUM

Pince à chalazion medium
Chalazion Klemme Medium

Overall Length : 90mm
Handle Length : 54 mm
21 mm diameter oval plate




CURETTE VARIOUS SIZES / Curettes à chalazion / Küretten verschieden Größen

Length : 110 mm
Handle length : 80 mm
CR01 Curette small
CR02 Curette medium,
(2.5 mm diameter chalazion curette)
CR03 Curette large,
(4 mm oval shaped chalazion curette)

CR01
CR02
CR03



CURETTES




CURETTE DOUBLE ENDED / Curette double / Kürette zweiseitig

Length : 130 mm
Handle length : 65 mm
Chalazion curette 2.5 mm diameter end
and oval end

CR05



CURETTES




CHALAZION CLAMP SMALL / Pince à chalazion 12 mm / Chalazion Klemme klein

Length : 115 mm
Handle length : 82 mm
12 mm diameter plate

CS01



CLAMPS



CHALAZION CLAMP LARGE / Pince à chalazion large, ovale / Chalazion Klemme groß

Length : 85 mm
Handle length : 53 mm
21 mm diameter oval plate

CS02



CLAMPS

