

CORNEAL & REFRACTIVE

Cornée & Réfractive
Korneal & Refraktiv


LASIK FLAP-LIFTING FORCEPS

Pince angulée mors long
Lasik Flap Hebepinzette

Overall length : 103 mm
Handle length : 78 mm
11 mm platform




SR14




REINSTEIN LENTICULE SEPARATOR

Séparateur Relex
Relex Häkchen

Length : 110 mm
Handle length : 90 mm
90° angled



SR70



DMEK BERNHEIM HOOK / Crochet DMEK de Bernheim / DMEK Häkchen nach Bernheim

Length : 120 mm
Handle length : 95 mm

SR109



CORNEAL EPITHELIUM SPATULA ROUND / Spatule pour épithélium cornéen, ronde / Epithelium Hornhautspatel, rund

Length : 110 mm
Handle length : 85 mm


SR100



CK LASIK MARKER / Marqueur LASIK avec croix de centrage / Hornhautmarker mit Fadenkreuz

Length : 100 mm
Handle length : 76 mm
Facilitate marking of the cornea
Round bodied

SR40



EPI PEELER / Scarificateur cornéen / Epi Peeler

Length : 110 mm
Handle length : 95 mm
Manipulate epithelium during surgery


SR51



EPI PEELER RECTANGLE SHAPED / Scarificateur, forme D / Epi Peeler - quadratische Form (LASIK)

Length : 110 mm
Handle length : 95 mm


SR61



EPI PEELER D SHAPED / Scarificateur Micro-Houe / Epi Peeler, gebogen

Length : 110 mm
Handle length : 95 mm


SR60



INTRA LASIK FLAP SEPARATOR / Elévateur volet cornéen, crochet pointu / Intra-Lasik Flap Separator

Length : 120 mm
Round bodied


SR34



KOCABA DMEK FORCEPS / Pince DMEK de Kocaba / Kocaba DMEK Pinzette

Length : 86 mm
Handle length : 56 mm


SR102



ALCOHOL WELL VARIOUS SIZES / Anneau marqueur / Lasek Alkohol Trichter

Length : 125 mm
Handle length : 92 mm
Liquid holding well for alcohol application
Round bodied
SR41 9 mm
SR42 8 mm
SR43 7 mm


SR41



LASIK CORNEAL MARKER / Marqueur cornéen en fourche / Lasik Hornhautmarker

Length : 120 mm
Handle length : 91 mm
Y shaped forked marker


SR21



LASIK FLAP ELEVATOR / Spatule Flap Lasik / LASIK Flap Heber

Length : 120 mm
Handle length : 95 mm
45° angled, 12 mm tip with 0.1 mm rounded end


SR31



LASIK FLAP FORCEPS / Pince Flap Lasik / Lasik Flap Pinzette

Length : 100 mm
Handle length : 77 mm
40° angled platform, circular 1 mm tips


SR71



LASIK FLAP SEPERATOR AND ELEVATOR DOUBLE END / Séparateur-élevateur LASIK, droit / Flap Heber gerade

Length : 120 mm
Handle length : 70 mm
45° and 65° ends
Double ended


SR103



LENTICULE SEPERATOR HOOK AND SPATULA DOUBLE END / Relex Crochet / Séparateur, Spatule / Relex Håkchen / Separator (Kombi-Instrument)

Length : 110 mm
Handle length : 70 mm
75° ends
Double ended


SR104



S MARKER / Marqueur DMEK «S» / S Marker

Length : 110 mm
Handle length : 90 mm
45° angle 2 mm marking end

SR101



TREPHINE VARIOUS SIZES / Trépan cornéen / Trepane verschiedene Größen

SR105 17 mm
SR106 18 mm
SR107 19 mm
SR110 16 mm


SR105



LASIK FEMTO SEPERATOR AND ELEVATOR DOUBLE ENDED / Séparateur-Élévateur Siebel / Lasik Femto Seperator und Elevator zweiseitig

Length : 115 mm
Handle length : 70.3 mm
Round bodied
Double ended


SR02



SINSKEY AND REVERSE SINSKEY HOOK / Manipulateur double Sinskey-Sinskey inversé / Sinskey Håkchen umgekel

Length : 120 mm
Handle length : 70.3 mm
Round bodied
Double ended

SR03



CORNEAL RING PUSHER / Manipulateur anneau intra-cornéen / Manipulator für Cornealen Ring

Length : 110 mm
Handle length : 91 mm
Round bodied

SR04



CORNEAL RING POCKETING HOOK / Crochet anneau intra-cornéen / Taschen Haken für Cornalen Ring

Length : 111 mm
Handle length : 91 mm
Round bodied


SR05



CORNEAL RING POCKET SPREADER / Spatule anneau intra-cornéen / Taschen Spreizer für Cornealen Ring

Overall Length : 108 mm
Handle Length : 90 mm


SR06



CORNEAL RING FORCEPS KERARING / Pince à anneau intra-cornéen - Keraring / Keraring Pinzette für Cornealen Ring

Overall Length : 103 mm
Handle Length : 77 mm


SR07



CORNEAL RING FORCEPS INTACS / Pince à anneau intra-cornéen - INTACS / Intacs Pinzette für Cornealen Ring

Length : 81 mm
Handle length : 53 mm


SR08



PATON SPATULA & SPOON DOUBLE ENDED / Spatule double de Paton / Paton Spatel und Löffel Zweiseitig


Length : 151 mm
Handle length : 71 mm
Round bodied
Double ended

SR108



LENTICULE EXTRACTION FORCEPS / Pince à extraction de lenticule / Lentikel Extraktion Pinzette

Overall Length : 115 mm
Handle Length : 81 mm



CORNEAL MARKER CAPSULORHEXIS GUIDE / Marqueur cornéen 6mm / Hornhautmarker für Kapsulhorexis

Length : 120 mm
Handle length : 92 mm
6 mm well
Round bodied

