




CORTEX REMOVAL




Irrigation





Aspiration




Male and female connectors are compatible with I/A and Phaco systems. Allows cortex removal through paracentesis or Phaco incision while maintaining the chamber with irrigation.



Aspiration



Hook dimensions 2.5mm long x 2.0mm wide



Hook dimensions 2.75mm long x 2.5mm wide

Bimanual I/A Cortex Removal System

10 per box

192-23

Irrigation 21GA x 1" (25mm)
Two .4mm ports with female adapter

Aspiration 23GA x 1" (25mm)
One .4mm port .5mm from tip
with male connector

Cortex Aspiration Cannula

10 per box

191-21

21GA x 1" (25mm) .4mm port with male adapter

191-23

23GA x 1" (25mm) .4mm port with male adapter

191-25

25GA x 1" (25mm) .4mm port with male adapter

Aspirating J-Hook Cannula

10 per box

104-25

25GA x 1" (25mm)

Aspirating J-Hook Cannula/ Capsule Polisher (Thornton)

10 per box

Open end port, roughened on the end of hook
for capsule polishing

107-25

25GA x 1" (25mm)


Made in USA

Qualitätsprodukte von 


salensa medizinischer Bedarf für das Auge

Fuggerstraße 1a • 88142 Wasserburg • Telefon 083 82 / 9430 225 • www.salensa.de • vertrieb@salensa.de


CORTEX REMOVAL




Hook dimensions 2.0mm long x 2.0mm wide



Hook dimensions 1.5mm long x 1.5mm wide



Hook dimensions 1.5mm long x 2.0mm wide



Aspirating J-Hook Cannula

10 per box

105-22L	22GA x 1" (25mm), Left
105-22R	22GA x 1" (25mm), Right
105-22S	22GA x 1" (25mm), Straight
105-25L	25GA x 1" (25mm), Left
105-25R	25GA x 1" (25mm), Right
105-25S	25GA x 1" (25mm), Straight
105-27L	27GA x 1" (25mm), Left
105-27R	27GA x 1" (25mm), Right
105-27S	27GA x 1" (25mm), Straight

Aspirating Cortex Extractor

10 per box

For aspirating cortex between the 10 and 2 o'clock positions

106-23S	23GA x 1" (25mm), Straight
106-25S	25GA x 1" (25mm), Straight
106-27L	27GA x 1" (25mm), Left
106-27R	27GA x 1" (25mm), Right
106-27S	27GA x 1" (25mm), Straight
106-30L	30GA x 1" (25mm), Left
106-30R	30GA x 1" (25mm), Right
106-30S	30GA x 1" (25mm), Straight

Aspirating C-Hook Extractor

10 per box

106-27SC	27GA
-----------------	------

Simcoe Cortex Extractor Aspiration Cannula 10 per box

Curved with .3mm inside port

108-23	23GA x 1" (25mm)
108-25	25GA x 1" (25mm)

Made in USA

Qualitätsprodukte von 