

GUARDED-SURGICAL KNIVES



Sub 2.0mm Phaco Slit Knife

6 per box

For bi-manual Phaco and micro-incision cataract procedures.

- EP-1014** 1.4mm, angled, full handle, bevel up
- EP-1015** 1.5mm, angled, full handle, bevel up
- EP-1016** 1.6mm, angled, full handle, bevel up
- EP-1017** 1.7mm, angled, full handle, bevel up
- EP-1018** 1.8mm, angled, full handle, bevel up
- EP-1020** 2.0mm, angled, full handle, bevel up



Clear Corneal Knife

6 per box

- EP-1122** 2.2 mm, angled, full handle, double bevel
- EP-1124** 2.4 mm, angled, full handle, double bevel
- EP-1127** 2.75 mm, angled, full handle, double bevel
- EP-1130** 3.0 mm, angled, full handle, double bevel
- EP-1132** 3.2 mm, angled, full handle, double bevel



Clear Corneal Knife, Parallel Sides

6 per box

- EP-1327** 2.75 mm, angled, full handle, double bevel
- EP-1330** 3.0 mm, angled, full handle, double bevel
- EP-1332** 3.2 mm, angled, full handle, double bevel



Trapezoid Clear Corneal / Implant Knife

6 per box

- EP-2832** 2.85 mm - 3.2 mm angled, full handle, bevel up
- EP-3032** 3.0 mm - 3.2 mm angled, full handle, bevel up

Dual-purpose design provides the ability to accomplish clear corneal and implant incisions with one knife.

GUARDED-SURGICAL KNIVES




Stab Knife

6 per box

For side port or scleral incision

EP-1500	15°, full handle
EP-2250	22.5°, full handle
EP-3000	30°, full handle
EP-4500	45°, full handle




Phaco Slit Knife

6 per box

For clear corneal or scleral incision

EP-5510	2.5 mm, angled, full handle, bevel up
EP-5568	2.65 mm, angled, full handle, bevel up
EP-4738	2.75 mm, angled, full handle, bevel up
EP-4740	2.85 mm, angled, full handle, bevel up
EP-4741	2.85 mm, straight, full handle
EP-5570	3.0 mm, angled, full handle, bevel up
EP-5571	3.0 mm, straight, full handle
EP-5572	3.2 mm, angled, full handle, bevel up
EP-5573	3.2 mm, straight, full handle
EP-5580	3.5 mm, angled, full handle, bevel up



Short Stop IOL Inserter Knife

6 per box

For IOL enlargement

EP-9336	3.5 mm, angled, full handle, bevel up
EP-9340	4.0 mm, angled, full handle, bevel up
EP-9350	5.2 mm, angled, full handle, bevel up




Crescent Knife

6 per box


For Scleral tunnel incision

EP-9402	2.0 mm, angled, full handle, bevel up
----------------	---------------------------------------


Instructions for using Guarded - Surgical Knives



Rotate the protective sleeve while gripping the end of the knife handle to release the "locking pin".



Slide the protective sleeve back while gripping the end of the knife handle.



Rotate the protective sleeve while gripping the end of the knife handle to secure the "locking pin" into open position.



Ready for use.

Reverse the above steps to secure the protective cover in the closed position prior to disposal into an approved sharps container.