

Ocular Irrigation Tube


20 per box

Extension tube with male/female luer connections provide fluid administration during ophthalmic procedures

169-30L-6 Irrigation tube - 6" (15.2cm)
(I.D. - .8mm, O.D. -1.6mm)

169-30L-10 Irrigation tube - 10" (25.4cm)
(I.D. - .8mm, O.D. -1.6mm)

169-30L-30 Irrigation tube - 30" (76.2cm)
(I.D. - 3mm, O.D. - 4mm)



Ocular Irrigation Set

3 per box

Features an adjustable roller clamp for control of fluid flow

Single

200-01 Tubing length 86"

Double

200-02 Tubing length 96"

Single w/ .20 micron filter

200-21 Tubing length 86"

Double w/ .20 micron filter

200-22 Tubing length 96"

Made in USA

Qualitätsprodukte von 