



GLAUCOMA



FLACS Cannula

5 per box


177-33 5mm beveled tip, slightly angled
Micro-gauge cannula for viscoelastic injection
during glaucoma surgery



Glaucoma Revision Pick

10 per box


60 Degree (Lederer)
177-27 27GA x 1/2" (13mm)
Angled 60°, 7mm bend to tip
Tip has a sharp side cutting edge
(superior approach)



Glaucoma Revision Pick

10 per box

30 Degree (Lederer)
178-27 27GA x 1/2" (13mm)
Angled 30°, 10mm bend to tip
Tip has a sharp side cutting edge
(temporal or nasal approach)



Anterior Chamber Cannula (Mendez)

10 per box

110-27TM 27GA x 10mm, round tip
Angled 45°, 3mm bend to tip to reform the anterior chamber
during glaucoma procedures