




# CORTEX REMOVAL



Irrigation



Aspiration




Male and female connectors are compatible with I/A and Phaco systems. Allows cortex removal through paracentesis or Phaco incision while maintaining the chamber with irrigation.


## Bimanual I/A Cortex Removal System **10 per box**

**192-23** Irrigation 21GA x 1" (25mm)  
Two .4mm ports with female adapter

Aspiration 23GA x 1" (25mm)  
One .4mm port .5mm from tip with male connector




Aspiration



## Cortex Aspiration Cannula **10 per box**



**191-21** 21GA x 1" (25mm) .4mm port with male adapter  
**191-23** 23GA x 1" (25mm) .4mm port with male adapter  
**191-25** 25GA x 1" (25mm) .4mm port with male adapter



Hook dimensions 2.5mm long x 2.0mm wide

## Aspirating J-Hook Cannula **10 per box**

**104-25** 25GA x 1" (25mm)





Hook dimensions 2.75mm long x 2.5mm wide

## Aspirating J-Hook Cannula/ Capsule Polisher (Thornton) **10 per box**

Open end port, roughened on the end of hook for capsule polishing

**107-25** 25GA x 1" (25mm)




Hook dimensions 2.0mm long x 2.0mm wide


## Aspirating J-Hook Cannula **10 per box**

**105-22L** 22GA x 1" (25mm), Left  
**105-22R** 22GA x 1" (25mm), Right  
**105-22S** 22GA x 1" (25mm), Straight  
**105-25L** 25GA x 1" (25mm), Left  
**105-25R** 25GA x 1" (25mm), Right  
**105-25S** 25GA x 1" (25mm), Straight  
**105-27L** 27GA x 1" (25mm), Left  
**105-27R** 27GA x 1" (25mm), Right  
**105-27S** 27GA x 1" (25mm), Straight


Made in USA

Qualitätsprodukte von 


# CORTEX REMOVAL



Hook dimensions 1.5mm long x 1.5mm wide



Hook dimensions 1.5mm long x 2.0mm wide



## Aspirating Cortex Extractor

10 per box

For aspirating cortex between the 10 and 2 o'clock positions

<b>106-23S</b>	23GA x 1" (25mm), Straight
<b>106-25S</b>	25GA x 1" (25mm), Straight
<b>106-27L</b>	27GA x 1" (25mm), Left
<b>106-27R</b>	27GA x 1" (25mm), Right
<b>106-27S</b>	27GA x 1" (25mm), Straight
<b>106-30L</b>	30GA x 1" (25mm), Left
<b>106-30R</b>	30GA x 1" (25mm), Right
<b>106-30S</b>	30GA x 1" (25mm), Straight

## Aspirating C-Hook Extractor

10 per box

<b>106-27SC</b>	27GA
-----------------	------

## Simcoe Cortex Extractor Aspiration Cannula 10 per box

Curved with .3mm inside port

<b>108-23</b>	23GA x 1" (25mm)
<b>108-25</b>	25GA x 1" (25mm)