



ANTERIOR CHAMBER MAINTENANCE



Anterior Chamber Maintainer Kit

3 per box


170-20T 20GA
(infusion cannula with self retaining serrated tip)
20GA x 3/16" (4.76mm)
cannula attached to 8" (20.3cm) flexible tubing
attached to an ocular irrigation set.



Anterior Chamber Maintainer

3 per box

170-20 20GA
(infusion cannula with self retaining serrated tip)
20GA x 3/16" (4.76mm) cannula
attached to 8" (20.3cm) flexible tubing



Pediatric A.C. Maintainer

3 per box

Used to maintain the chamber during pediatric retinal procedure (Mendez/Freeman)

171-20R 20GA
3.5mm, 90° angled, 20GA cannula
with 40° beveled tip connected to 11" (28cm)
flexible tube (.8mm I.D., 1.6mm O.D.)

Made in USA

Qualitätsprodukte von 