


# CAPSULE POLISHING




## Silicone "Squeegee" Tip Capsule Polisher (Engels)

10 per box

40 degree curve with extended silicone sleeve for polishing

**128-27** 27GA x 1" (25mm)




## Silicone "Squeegee" Tip Capsule Polisher (Terry)

10 per box

7mm bend to tip with flexible sleeve

**125-23** 23GA x 1" (25mm)

**125-27** 27GA x 1" (25mm)




## Silicone "Squeegee" Tip Cortex Extractor and Polisher (Buedding)

10 per box

Irrigating cortex extractor with flexible sleeve for the 12 o'clock position straight 23mm long

**107-27J** 27GA



## Capsule Polisher (Kratz)

10 per box

8mm bend, .30mm top port, rounded tip

**125-21** 21GA x 1" (25mm)

**125-24** 24GA x 1" (25mm)

Roughened end permits scraping and aspiration of cortical debris




## Capsule Polisher (Kratz)

10 per box


Roughened tip, angled 3mm

**125-22** 22GA x 1 1/2" (38mm)

Made in USA

Qualitätsprodukte von 

# CAPSULE POLISHING




## Capsule Polisher (Lu)

10 per box

8mm bend, .30mm top port, rounded tip

**124-23** 23GA x 1" (25mm)

Curved, roughened olive tip - open port




## Micro Loop Curette Polisher

10 per box

**125-25** 25GA x 1 1/4" (32mm)

Irrigating posterior capsule polisher  
curette loop is angled, posterior surface  
of loop is roughened, with open port



## Olive Tip Capsule Polisher

10 per box

10mm bend with roughened tip - open port

**126-23** 23GA x 1" (25mm)

**126-25** 25GA x 1" (25mm)