







CAPSULOTOMY



Irrigating Cystotome

10 per box

- 13mm bend to tip (formed)
120-01 25GA x 3/4" (19mm)
120-02 27GA x 3/4" (19mm)

Irrigating Cystotome/Push-Pull

10 per box

- 13mm bend to tip (formed)
120-03 25GA x 3/4" (19mm)

Irrigating Cystotome

10 per box

- Straight
120- 21GA x 1" (25mm)
120- 22GA x 1" (25mm)
120- 23GA x 1" (25mm)
120-25S 25GA x 5/8" (16mm)
120-27S 27GA x 1/2" (13mm)
120-30S 30GA x 3/4" (19mm)

Irrigating Cystotome

10 per box


- 15mm bend to tip (formed)
120-21F 21GA x 1" (25mm)
120-22F 22GA x 1" (25mm)
120-23F 23GA x 1" (25mm)
120-25F 25GA x 1" (25mm)

Irrigating Cystotome


10 per box

- 12mm bend to tip (formed)
120-22F-12 22GA x 3/4" (19mm)
120-23F-12 23GA x 3/4" (19mm)
120-25F-12 25GA x 5/8" (16mm)
120-27F-12 27GA x 1/2" (13mm)
 10mm bend to tip (formed)
120-30F-10 30GA x 1/2" (13mm)

Made in USA

Qualitätsprodukte von 

CAPSULOTOMY




Irrigating Cystotome

10 per box

8mm bend to tip (formed)

120-25F-8 25GA x 5/8" (16mm)

120-27F-8 27GA x 1/2" (13mm)




Irrigating Cystotome/Push-Pull (Blumenthal)

10 per box

Straight

121-25S 25GA x 5/8" (16mm)


121-27S 27GA x 1/2" (13mm)



8mm bend to tip (formed)

121-25F-8 25GA x 5/8" (16mm)


121-27F-8 27GA x 1/2" (13mm)



12mm bend to tip (formed)

121-25F-12 25GA x 5/8" (16mm)

121-27F-12 27GA x 1/2" (13mm)



Irrigating Cystotome/Reverse Bend


10 per box

Straight

122-23S 23GA x 3/4" (19mm)

122-25S 25GA x 5/8" (16mm)

122-27S 27GA x 1/2" (13mm)




12mm bend to tip (formed)

122-23RF 23GA x 3/4" (19mm)

122-25RF 25GA x 5/8" (16mm)

122-27RF 27GA x 1/2" (13mm)




Irrigating Cystotome (Joyce)

10 per box

Curved for intercapsular/endocapsular techniques

220-25 25GA x 5/8" (16mm)



Irrigating Cystotome (Joyce)


Curved for intercapsular/endocapsular techniques

220-27 27GA x 1/2" (13mm)

Made in USA

Qualitätsprodukte von 

CAPSULOTOMY




Capsulorhexis Suction Cystotome (Dahan) 10 per box

Curved design to create suction
and hold capsular flap

221-25* 25GA x 1" (25mm)

*Used with extension tubing #169-30L-10
and syringe part #300-03 or #300-05



Irrigating Cystotome/Angled 10 per box

Angled 45°, 10mm bend to tip

123-25 25GA x 5/8" (16mm)

123-27 27GA x 1/2" (13mm)