

# SURGICAL KNIVES



## Sub 2.0mm Phaco Slit Knife

6 per box

For bi-manual Phaco and micro-incision cataract procedures.

- EG-1014** 1.4mm, angled, full handle, bevel up
- EG-1015** 1.5mm, angled, full handle, bevel up
- EG-1016** 1.6mm, angled, full handle, bevel up
- EG-1017** 1.7mm, angled, full handle, bevel up
- EG-1018** 1.8mm, angled, full handle, bevel up
- EG-1020** 2.0mm, angled, full handle, bevel up



## Double Bevel Sideport Knife

6 per box

To create a self-sealing incision during bi-manual Phaco procedures.

- EG-5010** 1.0mm, angled, full handle, double bevel
- EG-5012** 1.2mm, angled, full handle, double bevel
- EG-5014** 1.4mm, angled, full handle, double bevel



## Spoon Knife

6 per box

For scleral tunnel incision.

- EG-6821** 3.0 mm, angled, full handle, double bevel



## MVR Knife

6 per box

For vitreoretinal infusion port or paracentesis incision.

- EG-5555** 19GA, full handle
- EG-5560** 20GA, full handle
- EG-5563** 23GA, full handle
- EG-5565** 25GA, full handle



## MVR Knife, Angled

6 per box

For vitreoretinal infusion port or paracentesis incision

- EG-5660** 20GA, full handle, angled
- EG-5663** 23GA, full handle, angled
- EG-5665** 25GA, full handle, angled

# SURGICAL KNIVES



## Clear Corneal Knife

6 per box

- EG-1122** 2.2 mm, angled, full handle, double bevel
- EG-1124** 2.4 mm, angled, full handle, double bevel
- EG-1127** 2.75 mm, angled, full handle, double bevel
- EG-1130** 3.0 mm, angled, full handle, double bevel
- EG-1132** 3.2 mm, angled, full handle, double bevel



## Clear Corneal Knife, Parallel Sides

6 per box

- EG-1327** 2.75 mm, angled, full handle, double bevel
- EG-1330** 3.0 mm, angled, full handle, double bevel
- EG-1332** 3.2 mm, angled, full handle, double bevel



## Trapezoid Clear Corneal / Implant Knife

6 per box

- EG-2832** 2.85 mm - 3.2 mm angled, full handle, bevel up
- EG-3032** 3.0 mm - 3.2 mm angled, full handle, bevel up

Dual-purpose design provides the ability to accomplish clear corneal and implant incisions with one knife.



## Stab Knife

6 per box

- EG-1500** 15°, full handle
- EG-2250** 22.5°, full handle
- EG-3000** 30°, full handle
- EG-4500** 45°, full handle



## Phaco Slit Knife

6 per box

For clear corneal or scleral incision

- EG-5510** 2.5 mm, angled, full handle, bevel up
- EG-5568** 2.65 mm, angled, full handle, bevel up
- EG-4738** 2.75 mm, angled, full handle, bevel up
- EG-4740** 2.85 mm, angled, full handle, bevel up
- EG-4741** 2.85 mm, straight, full handle
- EG-5570** 3.0 mm, angled, full handle, bevel up
- EG-5571** 3.0 mm, straight, full handle
- EG-5572** 3.2 mm, angled, full handle, bevel up
- EG-5573** 3.2 mm, straight, full handle
- EG-5580** 3.5 mm, angled, full handle, bevel up

Made in USA

Qualitätsprodukte von 

# SURGICAL KNIVES



## Short Stop IOL Inserter Knife

**6 per box**

For IOL enlargement

- EG-9336** 3.5 mm, angled, full handle, bevel up
- EG-9340** 4.0 mm, angled, full handle, bevel up
- EG-9350** 5.2 mm, angled, full handle, bevel up



## Crescent Knife

**6 per box**

For Scleral tunnel incision

- EG-9402** 2.0 mm, angled, full handle, bevel up

# GUARDED-SURGICAL KNIVES




## Sub 2.0mm Phaco Slit Knife

6 per box

For bi-manual Phaco and micro-incision cataract procedures.

- EP-1014** 1.4mm, angled, full handle, bevel up
- EP-1015** 1.5mm, angled, full handle, bevel up
- EP-1016** 1.6mm, angled, full handle, bevel up
- EP-1017** 1.7mm, angled, full handle, bevel up
- EP-1018** 1.8mm, angled, full handle, bevel up
- EP-1020** 2.0mm, angled, full handle, bevel up



## Clear Corneal Knife

6 per box

- EP-1122** 2.2 mm, angled, full handle, double bevel
- EP-1124** 2.4 mm, angled, full handle, double bevel
- EP-1127** 2.75 mm, angled, full handle, double bevel
- EP-1130** 3.0 mm, angled, full handle, double bevel
- EP-1132** 3.2 mm, angled, full handle, double bevel



## Clear Corneal Knife, Parallel Sides

6 per box

- EP-1327** 2.75 mm, angled, full handle, double bevel
- EP-1330** 3.0 mm, angled, full handle, double bevel
- EP-1332** 3.2 mm, angled, full handle, double bevel



## Trapezoid Clear Corneal / Implant Knife

6 per box

- EP-2832** 2.85 mm - 3.2 mm angled, full handle, bevel up
- EP-3032** 3.0 mm - 3.2 mm angled, full handle, bevel up

Dual-purpose design provides the ability to accomplish clear corneal and implant incisions with one knife.

# GUARDED-SURGICAL KNIVES




## Stab Knife

6 per box

For side port or scleral incision

- EP-1500** 15°, full handle
- EP-2250** 22.5°, full handle
- EP-3000** 30°, full handle
- EP-4500** 45°, full handle




## Phaco Slit Knife

6 per box

For clear corneal or scleral incision

- EP-5510** 2.5 mm, angled, full handle, bevel up
- EP-5568** 2.65 mm, angled, full handle, bevel up
- EP-4738** 2.75 mm, angled, full handle, bevel up
- EP-4740** 2.85 mm, angled, full handle, bevel up
- EP-4741** 2.85 mm, straight, full handle
- EP-5570** 3.0 mm, angled, full handle, bevel up
- EP-5571** 3.0 mm, straight, full handle
- EP-5572** 3.2 mm, angled, full handle, bevel up
- EP-5573** 3.2 mm, straight, full handle
- EP-5580** 3.5 mm, angled, full handle, bevel up



## Short Stop IOL Inserter Knife

6 per box

For IOL enlargement

- EP-9336** 3.5 mm, angled, full handle, bevel up
- EP-9340** 4.0 mm, angled, full handle, bevel up
- EP-9350** 5.2 mm, angled, full handle, bevel up




## Crescent Knife

6 per box

For Scleral tunnel incision

- EP-9402** 2.0 mm, angled, full handle, bevel up

# MICROSURGICAL BLADES



**Microsurgical Blade** **6 per box**

**EL-7511** 1.5mm, 15°, Orange



**EL-7512** 1.5mm, 15°, Orange



**Microsurgical Blade** **6 per box**

**EL-7513** 3.0mm, 15°, Blue



**EL-7514** 3.0mm, 15°, Blue



**Microsurgical Blade** **6 per box**

**EL-7515** 5.0mm, 15°, Green



**EL-7516** 5.0mm, 15°, Green



**Microsurgical Blade** **6 per box**

**EL-7522** 4.0mm, 22.5°, Grey



**EL-7523** 4.0mm, 22.5°, Grey



**Microsurgical Blade** **6 per box**

**EL-7530** 3.5mm, 30°, Black



**EL-7531** 3.5mm, 30°, Black




**Microsurgical Blade** **6 per box**

**EL-7545** 3.0mm, 45°, Tan



**EL-7546** 3.0mm, 45°, Tan

Made in USA

Qualitätsprodukte von 

# SURGICAL BLADES



## Needle Blade 6 per box

**EL-5910** 3.0mm



**EL-5920** 5.0mm



**EL-5930** 7.0mm



**EL-5940** 1.5mm



## Lamellar Blade 6 per box

**EL-6600** 60°



## Scleral Blade 6 per box

**EL-6823H** Full handle



## Scleral Blade 6 per box

**EL-5700**



## Miniature Edged Blade 20 per box

**EL-6200** Chisel tip



**EL-6400** Round tip



**EL-6500** Lance tip



**EL-6700** Curved tip




## Miniature Edged Blade 6 per box

**EL-6300** Pointed tip



**EL-6900** Sharp all around


Made in USA

Qualitätsprodukte von 


# Instructions for using Guarded - Surgical Knives



**Rotate the protective sleeve while gripping the end of the knife handle to release the "locking pin".**



**Slide the protective sleeve back while gripping the end of the knife handle.**



**Rotate the protective sleeve while gripping the end of the knife handle to secure the "locking pin" into open position.**



**Ready for use.**

**Reverse the above steps to secure the protective cover in the closed position prior to disposal into an approved sharps container.**