


LACRIMAL CANNULA




Lacrimal Cannula


10 per box

Open end

262-26 26GA x 1/2" (12mm) straight



260-26 26GA x 1" (25mm) curved



160-21 21GA x 1 1/2" (38mm) curved
160-25 25GA x 1 1/2" (38mm) curved
160-26 26GA x 1 1/2" (38mm) curved

Closed end with side port



161-23 23GA x 1 1/2" (38mm) curved, .3mm port
161-25 25GA x 1 1/2" (38mm) curved, .3mm port
163-21 21GA x 1 1/2" (38mm) straight, .3mm port
163-25 25GA x 1 1/2" (38mm) straight, .3mm port





Lacrimal Cannula (Straight)

10 per box

162-23 19GA x 1" (25mm) with 1/2" (13mm) 23GA tip

LACRIMAL CANNULA




Lacrimal Intubation Set

3 per box

180-23 Probe 23GA x 7" (18cm)
12" (30.5cm) flexible tubing
O.D. .025" (.60mm), I.D. .012" (.30mm)

180-27 Probe 27GA x 4 1/2" (11.5cm)
12" flexible tubing
O.D. .025" (.60mm), I.D. .012" (.30mm)



Lacrimal Intubation Set, (Crawford)

3 per box

180-23T Probe 23GA x 7" (18cm), olive tip
12" (30.5cm) flexible tubing
O.D. .025" (.60mm), I.D. .012" (.30mm)

180-27T Probe 27GA x 4 1/2" (11.5cm), olive tip
12" (30.5cm) flexible tubing
O.D. .025" (.60mm), I.D. .012" (.30mm)



Lacrimal Intubation Set, (Guibor)

3 per box

180-23C Probe 23GA x 7" (18cm), curved
12" (30.5cm) flexible tubing
O.D. .025" (.60mm), I.D. .012" (.30mm)