


LACRIMAL CANNULA




Lacrimal Cannula

10 per box

Open end

- 160-21** 21GA x 1 1/2" (38mm) curved
- 160-25** 25GA x 1 1/2" (38mm) curved
- 160-26** 26GA x 1 1/2" (38mm) curved




Closed end with side port

- 161-23** 23GA x 1 1/2" (38mm) curved, 0.3mm port
- 161-25** 25GA x 1 1/2" (38mm) curved, 0.3mm port
- 163-21** 21GA x 1 1/2" (38mm) straight, 0.3mm port
- 163-25** 25GA x 1 1/2" (38mm) straight, 0.3mm port

Qualitätsprodukte von



Made in USA



Lacrimal Cannula (Straight)

10 per box

- M9275** 26 G x 3/4" 0.45 x 20mm
Lacrimal cannula. Short straight

Qualitätsprodukte von

Sterimedix®

Made in UK



Lacrimal Cannula

10 per box

- JD90387** 26G 0.45x32mm, curved



- JD90388** 25G 0.5x29mm, curved



- JD90393** 26G disposable, 30 mm long, curved

Made in Italy

e.janach®

LACRIMAL CANNULA



Made in Italy

Lacrimal Dilator

10 per box

JD2622.1/10 Lacrimal dilator disposable, sterile

Lacrimal Intubation Set

5 per box

JD2622 Lacrimal duct cleaning disposable set, sterile
cannula: 26G curved cannula, 33 mm long
dilatore: Disposable lacrimal dilator, fine

Qualitätsprodukte von **e.janach**[®]