


# Sugi® DARC\*

# Sugi® Eyedrain



## Sterile Disposable Application Rings and Absorbent Strips

Systemic side effects may inadvertently occur when administering mydriatics for examining the back of the eye. The DARC (Drug Application Ring for the Cornea) allows precise application of the solution and its ring shape reduces to a minimum the risks and stress for the patient.

### Advantages and Benefits:

#### Sugi® DARC\*

- Precise application to the cornea reduces any systemic effect
- Uniform application
- Reduces the use of mydriatics
- Customised sizes
- Long-term effectiveness
- Highly biocompatible due to the mixture of cotton and regenerated cellulose
- No fibre residue on the eye

#### Sugi® Eyedrain

- Quickly drains fluids from the operation site
- Very gentle drainage
- Permanently clear view of the operation site
- High tensile strength
- Strips easily retained in the corner of the eye

### Ophthalmology

#### APPLICATIONS

##### Sugi® DARC\*

- Administering mydriatics for dilating the pupils
- Ophthalmic screening of infants

##### Sugi® Eyedrain

Draining fluids during

- Cataract extraction
- Phacoemulsification
- Vitrectomy

<b>Sugi® DARC*</b> 25 x 1 piece	Ø 14 x 8 mm	REF 30652
<b>Sugi® DARC*</b> 25 x 1 piece	Ø 10 x 5 mm	REF 30653
<b>Sugi® Eyedrain</b> 20 x 1 piece	155 x 4 x 1 mm	REF 30622

